

CORTEX™

The Future of Dental Implants



Product
Mini-Catalogue


Internal HEX Platform

Classix Prime Package									
Diameter Ø mm	Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.3	2.45	1.8	2.75			BCSC-1033	BCSC-1133	BCSC-1333	BCSC-1633
3.8	2.45	2.2	3.2		BCSC-0838	BCSC-1038	BCSC-1138	BCSC-1338	BCSC-1638
4.2	2.45	2.6	3.6	BCSC-0642	BCSC-0842	BCSC-1042	BCSC-1142	BCSC-1342	BCSC-1642
5.0	2.45	3.4	4.4	BCSC-0650	BCSC-0850	BCSC-1050	BCSC-1150	BCSC-1350	
6.0	2.45	4.4	5.4	BCSC-0660	BCSC-0860	BCSC-1060	BCSC-1160	BCSC-1360	

*The same prosthetic platform. C=Ø 3.65mm for all implant diameters

Classix Premium Package									
Diameter Ø mm	Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.3	2.45	1.8	2.75			PCSC-1033	PCSC-1133	PCSC-1333	PCSC-1633
3.8	2.45	2.2	3.2		PCSC-0838	PCSC-1038	PCSC-1138	PCSC-1338	PCSC-1638
4.2	2.45	2.6	3.6	PCSC-0642	PCSC-0842	PCSC-1042	PCSC-1142	PCSC-1342	PCSC-1642
5.0	2.45	3.4	4.4	PCSC-0650	PCSC-0850	PCSC-1050	PCSC-1150	PCSC-1350	
6.0	2.45	4.4	5.4	PCSC-0660	PCSC-0860	PCSC-1060	PCSC-1160	PCSC-1360	

*The same prosthetic platform. C=Ø 3.65mm for all implant diameters




Dynamix Prime Package									
Diameter Ø mm	Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.3	2.45	1.6	3.0			BCSD-1033	BCSD-1133	BCSD-1333	BCSD-1633
3.8	2.45	2.0	3.4		BCSD-0838	BCSD-1038	BCSD-1138	BCSD-1338	BCSD-1638
4.2	2.45	2.4	3.8	BCSD-0642	BCSD-0842	BCSD-1042	BCSD-1142	BCSD-1342	BCSD-1642
5.0	2.45	3.2	4.6	BCSD-0650	BCSD-0850	BCSD-1050	BCSD-1150	BCSD-1350	
6.0	2.45	4.2	5.6	BCSD-0660	BCSD-0860	BCSD-1060	BCSD-1160	BCSD-1360	

*The same prosthetic platform. C=Ø 3.65mm for all implant diameters

Dynamix Premium Package									
Diameter Ø mm	Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.3	2.45	1.6	3.0			PCSD-1033	PCSD-1133	PCSD-1333	PCSD-1633
3.8	2.45	2.0	3.4		PCSD-0838	PCSD-1038	PCSD-1138	PCSD-1338	PCSD-1638
4.2	2.45	2.4	3.8	PCSD-0642	PCSD-0842	PCSD-1042	PCSD-1142	PCSD-1342	PCSD-1642
5.0	2.45	3.2	4.6	PCSD-0650	PCSD-0850	PCSD-1050	PCSD-1150	PCSD-1350	
6.0	2.45	4.2	5.6	PCSD-0660	PCSD-0860	PCSD-1060	PCSD-1160	PCSD-1360	

*The same prosthetic platform. C=Ø 3.65mm for all implant diameters



Saturn Prime Package


Diameter Ø mm	Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.8	2.45	2.0	3.4	5.3	BST-0838	BST-1038	BST-1138	BST-1338	BST-1638
4.2	2.45	2.4	3.8	6	BST-0842	BST-1042	BST-1142	BST-1342	BST-1642

*The same prosthetic platform. C=Ø 3.65mm for all implant diameters

Saturn Premium Package


Diameter Ø mm	Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.8	2.45	2.0	3.4	5.3	PST-0838	PST-1038	PST-1138	PST-1338	PST-1638
4.2	2.45	2.4	3.8	6	PST-0842	PST-1042	PST-1142	PST-1342	PST-1642

*The same prosthetic platform. C=Ø 3.65mm for all implant diameters



1Piece Premium Package Smart

Diameter Ø mm	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.0	CP-1030	CP-1130	CP-1330	CP-1630
3.3	CP-1033	CP-1133	CP-1333	CP-1633



Conical Platform

Classix Conical Regular Platform - Prime Package


Diameter Ø [mm]	D.Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.3	2.1	2.0	3.2			BMCC-1033	BMCC-1133	BMCC-1333	BMCC-1633
3.8	2.1	2.5	3.7		BMCC-0838	BMCC-1038	BMCC-1138	BMCC-1338	BMCC-1638
4.2	2.1	2.7	4.0	BMCC-0642	BMCC-0842	BMCC-1042	BMCC-1142	BMCC-1342	BMCC-1642

*The same prosthetic platform. C=Ø 2.8mm for all implant diameters above

Classix Conical Wide Platform - Prime Package

Diameter Ø [mm]	D.Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
5.0	2.5	3.5	4.9	BMCC-0650	BMCC-0850	BMCC-1050	BMCC-1150	BMCC-1350	
6.0	2.5	4.4	5.9	BMCC-0660	BMCC-0860	BMCC-1060	BMCC-1160		

*The same prosthetic platform. C=Ø 3.8mm for all implant diameters above



Dynamix Conical Narrow Platform - Prime Package

Diameter Ø [mm]	D.Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.0	1.7	1.3	2.3			BMCD-1030	BMCD-1130	BMCD-1330	BMCD-1630

*The same prosthetic platform. C=Ø 2.3mm for all implant diameters above

Dynamix Conical Regular Platform - Prime Package


Diameter Ø [mm]	D.Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.3	2.1	1.6	2.6			BMCD-1033	BMCD-1133	BMCD-1333	BMCD-1633
3.8	2.1	1.9	2.9		BMCD-0838	BMCD-1038	BMCD-1138	BMCD-1338	BMCD-1638
4.2	2.1	2.2	3.2	BMCD-0642	BMCD-0842	BMCD-1042	BMCD-1142	BMCD-1342	BMCD-1642

*The same prosthetic platform. C=Ø 2.8mm for all implant diameters above

Dynamix Conical Wide Platform - Prime Package

Diameter Ø [mm]	D.Hex [mm]	A [mm]	B [mm]	L = 6 [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
5.0	2.5	2.8	3.8	BMCD-0650	BMCD-0850	BMCD-1050	BMCD-1150	BMCD-1350	
6.0	2.5	3.8	4.6	BMCD-0660	BMCD-0860	BMCD-1060	BMCD-1160		


*The same prosthetic platform. C=Ø 3.8mm for all implant diameters above























Magix Conical Regular Platform - Prime Package

Diameter Ø [mm]	D.Hex [mm]	A [mm]	B [mm]	L = 8 [mm]	L = 10 [mm]	L = 11.5 [mm]	L = 13 [mm]	L = 16 [mm]
3.3	2.1	1.6	2.6	BCDL-0833	BCDL-1033	BCDL-1133	BCDL-1333	BCDL-1633
3.8	2.1	1.9	2.9	BCDL-0838	BCDL-1038	BCDL-1138	BCDL-1338	BCDL-1638
4.2	2.1	2.2	3.2	BCDL-0842	BCDL-1042	BCDL-1142	BCDL-1342	BCDL-1642












*The same prosthetic platform. C=Ø 2.8mm for all implant diameters above
































	Cat.No	Product Description	Material
	CO-5000	Plastic Abutment Non Hex	POM
	CO-5020	Plastic Abutment with Hex	POM
	CO-5010	Wide Plastic Abutment with Hex	POM
	CO-5030	Wide Plastic Abutment Non Hex	POM
	CO-5100	Straight Multi-Unit abutment H1	Ti6Al4V-ELI
	CO-5200	Straight Multi-Unit abutment H2	Ti6Al4V-ELI
	CO-5300	Straight Multi-Unit abutment H3	Ti6Al4V-ELI
	CO-5400	Straight Multi-Unit abutment H4	Ti6Al4V-ELI
	CO-5500	Straight Multi-Unit abutment H5	Ti6Al4V-ELI
	CO-5118	Angulated Multi-Unit abutment 18Deg. H1	Ti6Al4V-ELI
	CO-5218	Angulated Multi-Unit abutment 18Deg. H2	Ti6Al4V-ELI
	CO-5318	Angulated Multi-Unit abutment 18Deg. H3	Ti6Al4V-ELI
	CO-5418	Angulated Multi-Unit abutment 18Deg. H4	Ti6Al4V-ELI
	CO-5518	Angulated Multi-Unit abutment 18Deg. H5	Ti6Al4V-ELI
	CO-5130	Angulated Multi-Unit abutment 30Deg. H1	Ti6Al4V-ELI
	CO-5230	Angulated Multi-Unit abutment 30Deg. H2	Ti6Al4V-ELI
	CO-5330	Angulated Multi-Unit abutment 30Deg. H3	Ti6Al4V-ELI
	CO-5430	Angulated Multi-Unit abutment 30Deg. H4	Ti6Al4V-ELI
	CO-5530	Angulated Multi-Unit abutment 30Deg. H5	Ti6Al4V-ELI
	CO-6020	POM (Polyacetal) for ball 2.5mm-Standard	POM
	CO-6021	POM (Polyacetal) for ball 2.5mm-Soft	POM
	CO-6022	POM (Polyacetal) for ball 2.5mm-Extra Soft	POM
	CO-6030	Metal Cap	Stainless Steel
	CO-6105	Aesthetic Connection Abutment L0.5mm	Ti6Al4V-ELI
	CO-6115	Aesthetic Connection Abutment L1.5mm	Ti6Al4V-ELI
	CO-6125	Aesthetic Connection Abutment L2.5mm	Ti6Al4V-ELI
	CO-6201	Ball Attachment 1mm	Ti6Al4V-ELI
	CO-6202	Ball Attachment 2mm	Ti6Al4V-ELI
	CO-6203	Ball Attachment 3mm	Ti6Al4V-ELI
	CO-6204	Ball Attachment 4mm	Ti6Al4V-ELI
	CO-6205	Ball Attachment 5mm	Ti6Al4V-ELI
	CO-6206	Ball Attachment 6mm	Ti6Al4V-ELI
	CO-6605	Plastic sleeve for Aesthetic Connection Abutment	POM
	CO-6606	Plastic sleeve for Aesthetic Connent Abut.Non Hex	POM



Cat.No	Product Description	Material
	CO-7102 Healing Cap 2mm	Ti6Al4V-ELI
	CO-7103 Healing Cap 3mm	Ti6Al4V-ELI
	CO-7104 Healing Cap 4mm	Ti6Al4V-ELI
	CO-7105 Healing Cap 5mm	Ti6Al4V-ELI
	CO-7106 Healing Cap 6mm	Ti6Al4V-ELI
	CO-7107 Healing Cap 7mm	Ti6Al4V-ELI
	CO-7102C Concave Hex. Healing Cap 2mm	Ti6Al4V-ELI
	CO-7103C Concave Hex. Healing Cap 3mm	Ti6Al4V-ELI
	CO-7104C Concave Hex. Healing Cap 4mm	Ti6Al4V-ELI
	CO-7105C Concave Hex. Healing Cap 5mm	Ti6Al4V-ELI
	CO-7106C Concave Hex. Healing Cap 6mm	Ti6Al4V-ELI
	CO-7503 Healing Cap 3mm Wide	Ti6Al4V-ELI
	CO-7504 Healing Cap 4mm Wide	Ti6Al4V-ELI
	CO-7505 Healing Cap 5mm Wide	Ti6Al4V-ELI
	CO-7506 Healing Cap 6mm Wide	Ti6Al4V-ELI
	CO-8000 Titanium Abutment L9mm	Ti6Al4V-ELI
	CO-8100 Titanium Abutment L13mm	Ti6Al4V-ELI
	CO-8000C Concave Straight Hex Abutment L9mm	Ti6Al4V-ELI
	CO-8100C Concave Straight Hex Abutment L13mm	Ti6Al4V-ELI
	CO-8028 Titanium connector abutment slim 5mm	Ti6Al4V-ELI
	CO-8029 Titanium connector abutment slim 7mm	Ti6Al4V-ELI
	CO-8031 Titanium connector abutment slim 9mm	Ti6Al4V-ELI
	CO-8036 Titanium Premium Connector Abutment	Ti6Al4V-ELI
	CO-8037 Titanium Connector Abutment 3mm collar	Ti6Al4V-ELI
	CO-8111 Titanium Abutment collar 1mm	Ti6Al4V-ELI
	CO-8112 Titanium Abutment collar 2mm	Ti6Al4V-ELI
	CO-8113 Titanium Abutment collar 3mm	Ti6Al4V-ELI
	CO-8114 Titanium Abutment collar 4mm	Ti6Al4V-ELI
	CO-8200 Titanium Abutment Ø 6.0mm L11mm	Ti6Al4V-ELI
	CO-8015 Angulated Titanium Abutment 15Deg. L9mm	Ti6Al4V-ELI
	CO-8115 Angulated Titanium abutment 15Deg. L12mm	Ti6Al4V-ELI
	CO-8025 Angulated Titanium Abutment 25Deg. L9mm	Ti6Al4V-ELI
	CO-8125 Angulated Titanium Abutment 25Deg. L12mm	Ti6Al4V-ELI









	Cat.No	Product Description	Material
	CO-8032	Angulated Titanium abutment slim 15Deg.	Ti6Al4V-ELI
	CO-8033	Angulated Titanium abutment slim 25Deg.	Ti6Al4V-ELI
	CO-8040	Implant Analog 3.75mm	Ti6Al4V-ELI
	CO-8047	Connector abutment analog	Stainless Steel
	CO-8049	Analog-Smart One Piece Square	Ti6Al4V-ELI
	CO-8050	Implant Transfer 9mm - Close Tray	Stainless Steel
	CO-8055	Implant Transfer 15mm - Close Tray	Stainless Steel
	CO-8060	Indirect Transfer 9mm open tray	Ti6Al4V-ELI
	CO-8065	Indirect Transfer 13mm open tray	Ti6Al4V-ELI
	CO-8060C	Concave Hex. Transfer 9mm open tray	Ti6Al4V-ELI
	CO-8065C	Concave Hex. Transfer 13mm open tray	Ti6Al4V-ELI
	CO-8070	Indirect narrow transfer for closed tray L9mm	Ti6Al4V-ELI
	CO-8075	Indirect narrow transfer for closed tray L13mm	Ti6Al4V-ELI
	CO-8071	Angulated Transfer Hexed Platform 15Deg. L9mm	Ti6Al4V-ELI
	CO-8072	Angulated Transfer Hexed Platform 25Deg. L9mm	Ti6Al4V-ELI
	CO-8080	Indirect standard transfer for closed tray L9mm	Ti6Al4V-ELI
	CO-8085	Indirect standart transfer for closed tray L13mm	Ti6Al4V-ELI
	CO-8080C	Concave Hex. Transfer for closed tray Platform 4.5mm L9mm	Ti6Al4V-ELI
	CO-8085C	Concave Hex. Transfer for closed tray Platform 4.5mm L13mm	Ti6Al4V-ELI
	CO-8102	Healing Cap L2mm Narrow	Ti6Al4V-ELI
	CO-8103	Healing Cap L3mm Narrow	Ti6Al4V-ELI
	CO-8104	Healing Cap L4mm Narrow	Ti6Al4V-ELI
	CO-8105	Healing Cap L5mm Narrow	Ti6Al4V-ELI
	CO-8106	Healing Cap L6mm Narrow	Ti6Al4V-ELI
	CO-8107	Healing Cap L7mm Narrow	Ti6Al4V-ELI
	CO-9001	Anatomic Titanium Abutment L1mm	Ti6Al4V-ELI
	CO-9002	Anatomic Titanium Abutment L2mm	Ti6Al4V-ELI
	CO-9003	Anatomic Titanium Abutment L3mm	Ti6Al4V-ELI
	CO-9004	Anatomic Titanium Abutment L4mm	Ti6Al4V-ELI
	CO-9001C	Anatomic Concave Straight Hex. Abutment L1mm	Ti6Al4V-ELI
	CO-9002C	Anatomic Concave Straight Hex. Abutment L2mm	Ti6Al4V-ELI
	CO-9003C	Anatomic Concave Straight Hex. Abutment L3mm	Ti6Al4V-ELI
	CO-9115	Anatomic Angulated Titanium Abutment 15Deg. L1mm	Ti6Al4V-ELI
	CO-9215	Anatomic Angulated Titanium Abutment 15Deg. L2mm	Ti6Al4V-ELI
	CO-9315	Anatomic Angulated Titanium Abutment 15Deg. L3mm	Ti6Al4V-ELI
	CO-9415	Anatomic Angulated Titanium Abutment 15Deg. L4mm	Ti6Al4V-ELI































	Cat.No	Product Description	Material
	CO-9115C	Anatomic Concave Angulated Hex. Abutment 15Deg. L1mm	Ti6Al4V-ELI
	CO-9215C	Anatomic Concave Angulated Hex. Abutment 15Deg. L2mm	Ti6Al4V-ELI
	CO-9315C	Anatomic Concave Angulated Hex. Abutment 15Deg. L3mm	Ti6Al4V-ELI
	CO-9125	Anatomic Angulated Titanium Abutment 25Deg. L1mm	Ti6Al4V-ELI
	CO-9225	Anatomic Angulated Titanium Abutment 25Deg. L2mm	Ti6Al4V-ELI
	CO-9325	Anatomic Angulated Titanium Abutment 25Deg. L3mm	Ti6Al4V-ELI
	CO-9425	Anatomic Angulated Titanium Abutment 25Deg. L4mm	Ti6Al4V-ELI
	CO-9125C	Anatomic Concave Angulated Hex. Abutment 25Deg. L1mm	Ti6Al4V-ELI
	CO-9225C	Anatomic Concave Angulated Hex. Abutment 25Deg. L2mm	Ti6Al4V-ELI
	CO-9325C	Anatomic Concave Angulated Hex. Abutment 25Deg. L3mm	Ti6Al4V-ELI
	CO-9700	Titanium Abutment Hex Dia 4.5mm	Ti6Al4V-ELI
	CO-9710	Titanium Abutment Without Hex Dia 4.5mm	Ti6Al4V-ELI
	CO-9701	Chrome-Cobalt Abutment Hex Dia 4.5mm	Cobalt Chromium
	CO-9711	Chrome-Cobalt Abutment Without Hex Dia 4.5mm	Cobalt Chromium
	CO-9701C	Chrome-Cobalt Concave Abutment Hex-Dia 4.5mm	Cobalt Chromium
	CO-9711C	Chrome-Cobalt Concave Abutment Without Hex-Dia 4.5mm	Cobalt Chromium
	CO-AA001	Anatomic Abutment for casting L1mm	POM
	CO-AA002	Anatomic Abutment for casting L2mm	POM
	CO-AA003	Anatomic Abutment for casting L3mm	POM
	CO-AA001C	Anatomic Concave Straight Castable Hex. Abutment L1mm	POM
	CO-AA002C	Anatomic Concave Straight Castable Hex. Abutment L2mm	POM
	CO-AA003C	Anatomic Concave Straight Castable Hex. Abutment L3mm	POM
	CO-AA151	Angulated abutment for casting 15Deg. L1mm	POM
	CO-AA152	Angulated abutment for casting 15Deg. L2mm	POM
	CO-AA153	Angulated abutment for casting 15Deg. L3mm	POM
	CO-AA151C	Anatomic Concave Angulated Castable Hex. Abutment 15Deg. L1mm	POM
	CO-AA152C	Anatomic Concave Angulated Castable Hex. Abutment 15Deg. L2mm	POM
	CO-AA153C	Anatomic Concave Angulated Castable Hex. Abutment 15Deg. L3mm	POM
	CO-AA251	Angulated abutment for casting 25Deg. L1mm	POM
	CO-AA252	Angulated abutment for casting 25Deg. L2mm	POM
	CO-AA253	Angulated abutment for casting 25Deg. L3mm	POM
	CO-AA251C	Anatomic Concave Angulated Castable Hex. Abutment 25Deg. L1mm	POM
	CO-AA252C	Anatomic Concave Angulated Castable Hex. Abutment 25Deg. L2mm	POM
	CO-AA253C	Anatomic Concave Angulated Castable Hex. Abutment 25Deg. L3mm	POM









	Cat.No	Product Description	Material
	CO-MUS65	Screw for angulated multi unit	Ti6Al4V-ELI
	CO-PA003	Temporary Peek Abutment L3mm	PEEK
	CO-PA153	Angulated Temporary Peek Abutment 15Deg. L3mm	PEEK
	CO-PA253	Angulated Temporary Peek Abutment 25Deg. L3mm	PEEK
	CO-PA003C	Anatomic Concave Straight Temporary Hex. Abutment L3mm	PEEK
	CO-PA153C	Anatomic Concave Angulated Temporary Hex. Abutment 15Deg. L3mm	PEEK
	CO-PA253C	Anatomic Concave Angulated Temporary Hex. Abutment 25Deg. L3mm	PEEK
	CO-S030	Abutment screw 7.8mm	Ti6Al4V-ELI
	CO-S060	Screw For Esthetic Connection Abutment L10.7mm	Ti6Al4V-ELI
	CO-S061	Screw For Esthetic Connection Abutment L11.3mm	Ti6Al4V-ELI
	CO-S062	Screw For Esthetic Connection Abutment L12.3mm	Ti6Al4V-ELI
	CO-SBMU	Scan Body for MU	Ti6Al4V-ELI
	CO-SBZZ	Scan Body INT/EXT- Hex	Ti6Al4V-ELI
	CO-SPCRC	ScanPost for CEREC HEX platform	Ti6Al4V-ELI
	CO-TBCRC	Ti-Base for CEREC Hex platform lock	Ti6Al4V-ELI
	CO-TBCRCR	Ti-Base for CEREC Hex platform no lock	Ti6Al4V-ELI
	CO-TBH1	T-base Hex. - H1mm L4mm lock	Ti6Al4V-ELI
	CO-TBH2	T-base Hex. - H1mm L6mm lock	Ti6Al4V-ELI
	CO-TBH3	T-base Hex. - H3mm L4mm lock	Ti6Al4V-ELI
	CO-TBH4	T-base Hex. - H3mm L6mm lock	Ti6Al4V-ELI
	CO-TBH1C	Concave T-base Hex. - H1mm L4mm lock	Ti6Al4V-ELI
	CO-TBH2C	Concave T-base Hex. - H1mm L6mm lock	Ti6Al4V-ELI
	CO-TBH3C	Concave T-base Hex. - H3mm L4mm lock	Ti6Al4V-ELI
	CO-TBH4C	Concave T-base Hex. - H3mm L6mm lock	Ti6Al4V-ELI
	CO-TBR1	Ti-Base Hex. -H1 L4. no lock	Ti6Al4V-ELI
	CO-TBR2	Ti-Base Hex. -H1 L6. no lock	Ti6Al4V-ELI
	CO-TBR3	Ti-Base Hex. -H3 L4. no lock	Ti6Al4V-ELI
	CO-TBR4	Ti-Base Hex. -H3 L6. no lock	Ti6Al4V-ELI
	CO-TBR1C	Concave T-base Hex. - H1mm L4mm no lock	Ti6Al4V-ELI
	CO-TBR2C	Concave T-base Hex. - H1mm L6mm no lock	Ti6Al4V-ELI
	CO-TBR3C	Concave T-base Hex. - H3mm L4mm no lock	Ti6Al4V-ELI
	CO-TBR4C	Concave T-base Hex. - H3mm L6mm no lock	Ti6Al4V-ELI
	R9000-R	Pre Mill Block - Hex	Ti6Al4V-ELI





















Cat.No	Product Description	Material
	MPN-5100 Straight Multi-Unit abutment H1 MPN	Ti6Al4V-ELI
	MPN-5200 Straight Multi-Unit abutment H2 MPN	Ti6Al4V-ELI
	MPN-5300 Straight Multi-Unit abutment H3 MPN	Ti6Al4V-ELI
	MPN-5400 Straight Multi-Unit abutment H4 MPN	Ti6Al4V-ELI
	MPN-5500 Straight Multi-Unit abutment H5 MPN	Ti6Al4V-ELI
	MPN-5118 Angulated Multi-Unit abutment 18Deg. H1 MPN	Ti6Al4V-ELI
	MPN-5218 Angulated Multi-Unit abutment 18Deg. H2 MPN	Ti6Al4V-ELI
	MPN-5318 Angulated Multi-Unit abutment 18Deg. H3 MPN	Ti6Al4V-ELI
	MPN-5418 Angulated Multi-Unit abutment 18Deg. H4 MPN	Ti6Al4V-ELI
	MPN-5518 Angulated Multi-Unit abutment 18Deg. H5 MPN	Ti6Al4V-ELI
	MPN-5130 Angulated Multi-Unit abutment 30Deg. H1 MPN	Ti6Al4V-ELI
	MPN-5230 Angulated Multi-Unit abutment 30Deg. H2 MPN	Ti6Al4V-ELI
	MPN-5330 Angulated Multi-Unit abutment 30Deg. H3 MPN	Ti6Al4V-ELI
	MPN-5430 Angulated Multi-Unit abutment 30Deg. H4 MPN	Ti6Al4V-ELI
	MPN-5530 Angulated Multi-Unit abutment 30Deg. H5 MPN	Ti6Al4V-ELI
	MPN-6201 Ball Attachment 1mm MPN	Ti6Al4V-ELI
	MPN-6202 Ball Attachment 2mm MPN	Ti6Al4V-ELI
	MPN-6203 Ball Attachment 3mm MPN	Ti6Al4V-ELI
	MPN-6204 Ball Attachment 4mm MPN	Ti6Al4V-ELI
	MPN-6205 Ball Attachment 5mm MPN	Ti6Al4V-ELI
	MPN-8000 Titanium Abutment L9mm MPN	Ti6Al4V-ELI
	MPN-8036 Premium Straight Abutment MPN	Ti6Al4V-ELI
	MPN-8040 Implant Analog MPN	Ti6Al4V-ELI
	MPN-8060 Indirect Transfer Open Tray MPN	Ti6Al4V-ELI
	MPN-8080 Impression transfer for close tray L9mm MPN	Ti6Al4V-ELI
	MPN-8085 Impression transfer for close tray L12mm MPN	Ti6Al4V-ELI
	MPN-8102 Healing Cap L2mm MPN	Ti6Al4V-ELI
	MPN-8103 Healing Cap L3mm MPN	Ti6Al4V-ELI
	MPN-8104 Healing Cap L4mm MPN	Ti6Al4V-ELI
	MPN-8105 Healing Cap L5mm MPN	Ti6Al4V-ELI
	MPN-8106 Healing Cap L6mm MPN	Ti6Al4V-ELI
	MPN-8107 Healing Cap L7mm MPN	Ti6Al4V-ELI
	MPN-8115 Angulated Titanium Abutment 15Deg. MPN	Ti6Al4V-ELI






















	Cat.No	Product Description	Material
	MPN-8125	Angulated Titanium Abutment 25Deg. MPN	Ti6Al4V-ELI
	MPN-9001	Anatomic Titanium Abutment L1mm MPN	Ti6Al4V-ELI
	MPN-9002	Anatomic Titanium Abutment L2mm MPN	Ti6Al4V-ELI
	MPN-9003	Anatomic Titanium Abutment L3mm MPN	Ti6Al4V-ELI
	MPN-9004	Anatomic Titanium Abutment L4mm MPN	Ti6Al4V-ELI
	MPN-9115	Anatomic Ang. Titanium Abutment 15Deg. L1mm MPN	Ti6Al4V-ELI
	MPN-9215	Anatomic Ang. Titanium Abutment 15Deg. L2mm MPN	Ti6Al4V-ELI
	MPN-9315	Anatomic Ang. Titanium Abutment 15Deg. L3mm MPN	Ti6Al4V-ELI
	MPN-9415	Anatomic Ang. Titanium Abutment 15Deg. L4mm MPN	Ti6Al4V-ELI
	MPN-9225	Anatomic Ang Titanium Abutment 25Deg. L2mm MPN	Ti6Al4V-ELI
	MPN-9325	Anatomic Ang. Titanium Abutment 25Deg. L3mm MPN	Ti6Al4V-ELI
	MPN-9425	Anatomic Ang. Titanium Abutment 25Deg. L4mm MPN	Ti6Al4V-ELI
	MPN-9700	Titanium Castable Abutment Hex MPN	Ti6Al4V-ELI
	MPN-9710	Titanium Castable Abutment Non Hex MPN	Ti6Al4V-ELI
	MPN-9701	CrCo Castable Abutment Hex MPN	Cobalt Chromium
	MPN-9711	CrCo Castable Abutment Non Hex MPN	Cobalt Chromium
	MPN-MUS65	Screw for Narrow Platform Multi-Unit Abutment	Ti6Al4V-ELI
	MPN-PA00	Straight Temporary Abutment L2mm MPN	PEEK
	MPN-PA15	Angulated Temporary Abutment 15Deg. L2mm MPN	PEEK
	MPN-PA25	Angulated Temporary Abutment 25Deg. L2mm MPN	PEEK
	MPN-S030	Abutment Screw MPN	Ti6Al4V-ELI
	MPN-SBZZ	Scan Body INT/EXT- Narrow Platform	Ti6Al4V-ELI
	MPN-SPCRC	ScanPost for CEREC Narrow platform	Ti6Al4V-ELI
	MPN-TBCRC	Ti-Base for CEREC Narrow platform lock	Ti6Al4V-ELI
	MPN-TBCRCR	Ti-Base for CEREC Narrow platform no lock	Ti6Al4V-ELI
	MPN-TBH1	Ti-Base H1 - Narrow Platform lock	Ti6Al4V-ELI
	MPN-TBH2	Ti-Base H2 - Narrow Platform lock	Ti6Al4V-ELI
	MPN-TBH3	Ti-Base H3 - Narrow Platform lock	Ti6Al4V-ELI
	MPN-TBR1	Ti-Base H1 - Narrow Platform no lock	Ti6Al4V-ELI
	MPN-TBR2	Ti-Base H2 - Narrow Platform no lock	Ti6Al4V-ELI
	MPN-TBR3	Ti-Base H3 - Narrow Platform no lock	Ti6Al4V-ELI
	S9030-R	Pre Mill Block - Narrow Platform	Ti6Al4V-ELI
	MPR-0055	Low Profile Abutment 5.5mm MPR	Ti6Al4V-ELI
	MPR-0070	Low Profile Abutment 7mm MPR	Ti6Al4V-ELI















Cat.No	Product Description	Material
	MPR-5100 Straight Multi-Unit abutment H1 MPR	Ti6Al4V-ELI
	MPR-5200 Straight Multi-Unit abutment H2 MPR	Ti6Al4V-ELI
	MPR-5300 Straight Multi-Unit abutment H3 MPR	Ti6Al4V-ELI
	MPR-5400 Straight Multi-Unit abutment H4 MPR	Ti6Al4V-ELI
	MPR-5500 Straight Multi-Unit abutment H5 MPR	Ti6Al4V-ELI
	MPR-5118 Angulated Multi-Unit abutment 18Deg. H1 MPR	Ti6Al4V-ELI
	MPR-5218 Angulated Multi-Unit abutment 18Deg. H2 MPR	Ti6Al4V-ELI
	MPR-5318 Angulated Multi-Unit abutment 18Deg. H3 MPR	Ti6Al4V-ELI
	MPR-5418 Angulated Multi-Unit abutment 18Deg. H4 MPR	Ti6Al4V-ELI
	MPR-5518 Angulated Multi-Unit abutment 18Deg. H5 MPR	Ti6Al4V-ELI
	MPR-5130 Angulated Multi-Unit abutment 30Deg. H1 MPR	Ti6Al4V-ELI
	MPR-5230 Angulated Multi-Unit abutment 30Deg. H2 MPR	Ti6Al4V-ELI
	MPR-5330 Angulated Multi-Unit abutment 30Deg. H3 MPR	Ti6Al4V-ELI
	MPR-5430 Angulated Multi-Unit abutment 30Deg. H4 MPR	Ti6Al4V-ELI
	MPR-5530 Angulated Multi-Unit abutment 30Deg. H5 MPR	Ti6Al4V-ELI
	MPR-6201 Ball Attachment 1mm MPR	Ti6Al4V-ELI
	MPR-6202 Ball Attachment 2mm MPR	Ti6Al4V-ELI
	MPR-6203 Ball Attachment 3mm MPR	Ti6Al4V-ELI
	MPR-6204 Ball Attachment 4mm MPR	Ti6Al4V-ELI
	MPR-6205 Ball Attachment 5mm MPR	Ti6Al4V-ELI
	MPR-6206 Ball Attachment 6mm MPR	Ti6Al4V-ELI
	MPR-6207 Ball Attachment 7mm MPR	Ti6Al4V-ELI
	MPR-7102 Healing Cap 2mm MPR	Ti6Al4V-ELI
	MPR-7103 Healing Cap 3mm MPR	Ti6Al4V-ELI
	MPR-7104 Healing Cap 4mm MPR	Ti6Al4V-ELI
	MPR-7105 Healing Cap 5mm MPR	Ti6Al4V-ELI
	MPR-7106 Healing Cap 6mm MPR	Ti6Al4V-ELI
	MPR-7107 Healing Cap 7mm MPR	Ti6Al4V-ELI
	MPR-7503 Healing Cap Dia. 6mm L3mm MPR	Ti6Al4V-ELI
	MPR-7505 Healing Cap Dia. 6mm L5mm MPR	Ti6Al4V-ELI
	MPR-7507 Healing Cap Dia. 6mm L7mm MPR	Ti6Al4V-ELI
	MPR-8000 Titanium Abutment L9mm MPR	Ti6Al4V-ELI
	MPR-8036 Premium Straight Abutment MPR	Ti6Al4V-ELI



































	Cat.No	Product Description	Material
	MPR-8040	Implant Analog MPR	Ti6Al4V-ELI
	MPR-8060	Indirect Transfer Open Tray MPR	Ti6Al4V-ELI
	MPR-8080	Impression transfer for close tray L9mm MPR	Ti6Al4V-ELI
	MPR-8085	Impression transfer for close tray L12mm MPR	Ti6Al4V-ELI
	MPR-8115	Angulated Titanium Abutment 15Deg. MPR	Ti6Al4V-ELI
	MPR-8125	Angulated Titanium Abutment 25Deg. MPR	Ti6Al4V-ELI
	MPR-9001	Anatomic Titanium Abutment L1mm MPR	Ti6Al4V-ELI
	MPR-9002	Anatomic Titanium Abutment L2mm MPR	Ti6Al4V-ELI
	MPR-9003	Anatomic Titanium Abutment L3mm MPR	Ti6Al4V-ELI
	MPR-9004	Anatomic Titanium Abutment L4mm MPR	Ti6Al4V-ELI
	MPR-9115	Anatomic Ang. Titanium Abutment 15Deg. L1mm MPR	Ti6Al4V-ELI
	MPR-9215	Anatomic Ang. Titanium Abutment 15Deg. L2mm MPR	Ti6Al4V-ELI
	MPR-9315	Anatomic Ang. Titanium Abutment 15Deg. L3mm MPR	Ti6Al4V-ELI
	MPR-9415	Anatomic Ang. Titanium Abutment 15Deg. L4mm MPR	Ti6Al4V-ELI
	MPR-9225	Anatomic Ang. Titanium Abutment 25Deg. L2mm MPR	Ti6Al4V-ELI
	MPR-9325	Anatomic Ang. Titanium Abutment 25Deg. L3mm MPR	Ti6Al4V-ELI
	MPR-9425	Anatomic Ang. Titanium Abutment 25Deg. L4mm MPR	Ti6Al4V-ELI
	MPR-9700	Titanium Castable Abutment Hex MPR	Ti6Al4V-ELI
	MPR-9710	Titanium Castable Abutment Non Hex MPR	Ti6Al4V-ELI
	MPR-9701	CrCo Castable Abutment Hex MPR	Cobalt Chromium
	MPR-9711	CrCo Castable Abutment Non Hex MPR	Cobalt Chromium
	MPR-MUS65	Screw for Regular Platform Multi-Unit Abutment	Ti6Al4V-ELI
	MPR-PA00	Straight Temporary Abutment L2mm MPR	PEEK
	MPR-PA15	Angulated Temporary Abutment 15Deg. L2mm MPR	PEEK
	MPR-PA25	Angulated Temporary Abutment 25Deg. L2mm MPR	PEEK
	MPR-S030	Abutment Screw MPR	Ti6Al4V-ELI
	MPR-SBZZ	Scan Body INT/EXT- Regular Platform	Ti6Al4V-ELI
	MPR-SPCRC	ScanPost for CEREC Regular platform	Ti6Al4V-ELI
	MPR-TBCRC	Ti-Base for CEREC Regular platform lock	Ti6Al4V-ELI
	MPR-TBCRCR	Ti-Base for CEREC Regular platform no lock	Ti6Al4V-ELI

























Cat.No	Product Description	Material		
	MPR-TBH1	Ti-Base H1 Regular Platform lock	Ti6Al4V-ELI	
	MPR-TBH2	Ti-Base H2 Regular Platform lock	Ti6Al4V-ELI	
	MPR-TBH3	Ti-Base H3 Regular Platform lock	Ti6Al4V-ELI	
	MPR-TBH4	Ti-Base H4 Regular Platform lock	Ti6Al4V-ELI	
	MPR-TBR1	Ti-Base H1 Regular Platform no lock	Ti6Al4V-ELI	
	MPR-TBR2	Ti-Base H2 Regular Platform no lock	Ti6Al4V-ELI	
	MPR-TBR3	Ti-Base H3 Regular Platform no lock	Ti6Al4V-ELI	
	MPR-TBR4	Ti-Base H4 Regular Platform no lock	Ti6Al4V-ELI	
	S9000-R	Pre Mill Block - Regular Platform	Ti6Al4V-ELI	
	MPW-5100	Straight Multi-Unit abutment H1 MPW	Ti6Al4V-ELI	
	MPW-5200	Straight Multi-Unit abutment H2 MPW	Ti6Al4V-ELI	
	MPW-5300	Straight Multi-Unit abutment H3 MPW	Ti6Al4V-ELI	
	MPW-5400	Straight Multi-Unit abutment H4 MPW	Ti6Al4V-ELI	
	MPW-5500	Straight Multi-Unit abutment H5 MPW	Ti6Al4V-ELI	
		MPW-5118	Angulated Multi-Unit abutment 18Deg. H1 MPW	Ti6Al4V-ELI
		MPW-5218	Angulated Multi-Unit abutment 18Deg. H2 MPW	Ti6Al4V-ELI
		MPW-5318	Angulated Multi-Unit abutment 18Deg. H3 MPW	Ti6Al4V-ELI
		MPW-5418	Angulated Multi-Unit abutment 18Deg. H4 MPW	Ti6Al4V-ELI
		MPW-5518	Angulated Multi-Unit abutment 18Deg. H5 MPW	Ti6Al4V-ELI
	MPW-5130	Angulated Multi-Unit abutment 30Deg. H1 MPW	Ti6Al4V-ELI	
	MPW-5230	Angulated Multi-Unit abutment 30Deg. H2 MPW	Ti6Al4V-ELI	
	MPW-5330	Angulated Multi-Unit abutment 30Deg. H3 MPW	Ti6Al4V-ELI	
	MPW-5430	Angulated Multi-Unit abutment 30Deg. H4 MPW	Ti6Al4V-ELI	
	MPW-5530	Angulated Multi-Unit abutment 30Deg. H5 MPW	Ti6Al4V-ELI	
	MPW-6201	Ball Attachment 1mm MPW	Ti6Al4V-ELI	
	MPW-6202	Ball Attachment 2mm MPW	Ti6Al4V-ELI	
	MPW-6203	Ball Attachment 3mm MPW	Ti6Al4V-ELI	
	MPW-6204	Ball Attachment 4mm MPW	Ti6Al4V-ELI	
	MPW-6205	Ball Attachment 5mm MPW	Ti6Al4V-ELI	
	MPW-6206	Ball Attachment 6mm MPW	Ti6Al4V-ELI	
	MPW-7503	Healing Cap 3mm MPW	Ti6Al4V-ELI	
	MPW-7504	Healing Cap 4mm MPW	Ti6Al4V-ELI	
	MPW-7505	Healing Cap 5mm MPW	Ti6Al4V-ELI	
	MPW-7506	Healing Cap 6mm MPW	Ti6Al4V-ELI	

	Cat.No	Product Description	Material
	MPW-7603	Healing Cap Dia7mm L3mm MPW	Ti6Al4V-ELI
	MPW-7605	Healing Cap Dia7mm L5mm MPW	Ti6Al4V-ELI
	MPW-7607	Healing Cap Dia7mm L7mm MPW	Ti6Al4V-ELI
	MPW-8000	Titanium Abutment L9mm MPW	Ti6Al4V-ELI
	MPW-8036	Premium Straight Abutment MPW	Ti6Al4V-ELI
	MPW-8040	Implant Analog MPW	Ti6Al4V-ELI
	MPW-8060	Indirect Transfer Open tray MPW	Ti6Al4V-ELI
	MPW-8080	Impression transfer for close tray L9mm MPW	Ti6Al4V-ELI
	MPW-8085	Impression transfer for close tray 12mm MPW	Ti6Al4V-ELI
	MPW-8115	Angulated Titanium Abutment 15Deg. MPW	Ti6Al4V-ELI
	MPW-8125	Angulated Titanium Abutment 25Deg. MPW	Ti6Al4V-ELI
	MPW-9001	Anatomic Titanium Abutment L1mm MPW	Ti6Al4V-ELI
	MPW-9002	Anatomic Titanium Abutment L2mm MPW	Ti6Al4V-ELI
	MPW-9003	Anatomic Titanium Abutment L3mm MPW	Ti6Al4V-ELI
	MPW-9004	Anatomic Titanium Abutment L4mm MPW	Ti6Al4V-ELI
	MPW-9115	Anatomic Ang. Titanium Abutment 15Deg. L1mm MPW	Ti6Al4V-ELI
	MPW-9215	Anatomic Ang. Titanium Abutment 15Deg. L2mm MPW	Ti6Al4V-ELI
	MPW-9315	Anatomic Ang. Titanium Abutment 15Deg. L3mm MPW	Ti6Al4V-ELI
	MPW-9415	Anatomic Ang. Titanium Abutment 15Deg. L4mm MPW	Ti6Al4V-ELI
	MPW-9125	Anatomic Ang. Titanium Abutment 25Deg. L1mm MPW	Ti6Al4V-ELI
	MPW-9225	Anatomic Ang. Titanium Abutment 25Deg. L2mm MPW	Ti6Al4V-ELI
	MPW-9325	Anatomic Ang. Titanium Abutment 25Deg. L3mm MPW	Ti6Al4V-ELI
	MPW-9425	Anatomic Ang. Titanium Abutment 25Deg. L4mm MPW	Ti6Al4V-ELI
	MPW-9700	Titanium Castable Abutment Hex MPW	Ti6Al4V-ELI
	MPW-9710	Titanium Castable Abutment Non Hex MPW	Ti6Al4V-ELI
	MPW-9701	CrCo Castable Abutment Hex MPW	Cobalt Chromium
	MPW-9711	CrCo Castable Abutment Non Hex MPW	Cobalt Chromium
	MPW-MUS65	Screw for Wide Platform Multi-Unit Abutment	Ti6Al4V-ELI
	MPW-PA00	Straight Temporary Abutment L2mm MPW	PEEK
	MPW-PA15	Angulated Temporary Abutment 15Deg. L2mm MPW	PEEK
	MPW-PA25	Angulated Temporary Abutment 25Deg. L2mm MPW	PEEK
	MPW-S030	Abutment Screw MPW	Ti6Al4V-ELI
	MPW-SBZZ	Scan Body INT/EXT- Wide Platform	Ti6Al4V-ELI
	MPW-SPCRC	ScanPost for CEREC Wide platform	Ti6Al4V-ELI

	Cat.No	Product Description	Material
	MPW-TBCRC	Ti-Base for CEREC Wide platform lock	Ti6Al4V-ELI
	MPW-TBCRCR	Ti-Base for CEREC Wide platform no lock	Ti6Al4V-ELI
	MPW-TBH1	Ti-Base H1 Wide Platform lock	Ti6Al4V-ELI
	MPW-TBH2	Ti-Base H2 Wide Platform lock	Ti6Al4V-ELI
	MPW-TBH3	Ti-Base H3 Wide Platform lock	Ti6Al4V-ELI
	MPW-TBH4	Ti-Base H4 Wide Platform lock	Ti6Al4V-ELI
	MPW-TBR1	Ti-Base H1 Wide Platform no lock	Ti6Al4V-ELI
	MPW-TBR2	Ti-Base H2 Wide Platform no lock	Ti6Al4V-ELI
	MPW-TBR3	Ti-Base H3 Wide Platform no lock	Ti6Al4V-ELI
	MPW-TBR4	Ti-Base H4 Wide Platform no lock	Ti6Al4V-ELI
	S9020-R	Pre Mill Block - Wide Platform	Ti6Al4V-ELI
	CO-0620	POM (Polyacetal) Standard (Yellow)	POM
	CO-0621	POM (Polyacetal) Soft (White)	POM
	CO-0622	POM (Polyacetal) Hard (Red)	POM
	CO-0630	Metallic Attachment Housing	Ti6Al4V-ELI
	MU-5060	Transfer for Multi-Unit abutment open tray	Ti6Al4V-ELI
	MU-5080	Transfer for Multi-Unit abutment close tray	Ti6Al4V-ELI
	MU-50HC	Healing cap for Multi-Unit abutment	Ti6Al4V-ELI
	MU-50PL	Plastic sleeve for Multi-Unit abutment	POM
	MU-50SPL	Plastic sleeve for single unit abutment	POM
	MU-50TL	Titanium sleeve for Multi-Unit abutment	Ti6Al4V-ELI
	MU-50TS	Short Titanium sleeve for Multi-Unit abutment	Ti6Al4V-ELI
	MU-8041	Multi-Unit abutment analog	Ti6Al4V-ELI
	MU-STSSST	Short Ti-Base (sleeve) for MU abutment	Ti6Al4V-ELI
	MU-LTSST	Long Ti-Base (sleeve) for MU abutment	Ti6Al4V-ELI
	MU-S010	Short Screw for MU abutment	Ti6Al4V-ELI
	MU-S020	Screw for Multi-Unit abutment	Ti6Al4V-ELI
	MU-S060	Screw for Multi-Unit open tray transfer	Ti6Al4V-ELI
	MU-S080	Screw for Multi-Unit close tray transfer	Ti6Al4V-ELI
	CD-1001	Starter drill 1.50mm	Stainless Steel
	CD-5003	Pilot Drill Ø=2.00mm White	Stainless Steel

	Cat.No	Product Description	Material
	CD-5010	Drill Ø=2.8mm Yellow	Stainless Steel
	CD-5137	Drill for bone type D1 Dia.3.7mm for implant 3.8	Stainless Steel
	CD-5140	Drill for bone type D1 Dia.4.1mm for implant 4.2	Stainless Steel
	CD-5148	Drill for bone type D1 Dia.4.9mm for implant 5.0	Stainless Steel
	CD-C038	Conical Drill dia. 3.2mm for implant dia. 3.8mm	Stainless Steel
	CD-C042	Conical Drill dia. 3.7mm for implant dia. 4.2mm	Stainless Steel
	CD-C050	Conical Drill dia. 4.3mm for implant dia. 5.0mm	Stainless Steel
	CD-C060	Conical Drill dia. 5.4mm for implant dia. 6.0mm	Stainless Steel
	CD-GS1001	GS Start Drill	Stainless Steel
	CD-GS0822	GS Drill D2,0x8 mm DLC	Stainless Steel
	CD-GS1022	GS Drill D2,0x10 mm DLC	Stainless Steel
	CD-GS1122	GS Drill D2,0x11,5 mm DLC	Stainless Steel
	CD-GS1322	GS Drill D2,0x13 mm DLC	Stainless Steel
	CD-GS1622	GS Drill D2,0x16 mm DLC	Stainless Steel
	CD-GS0828	GS Drill D2,8x8 mm DLC	Stainless Steel
	CD-GS1028	GS Drill D2,8x10 mm DLC	Stainless Steel
	CD-GS1128	GS Drill D2,8x11,5 mm DLC	Stainless Steel
	CD-GS1328	GS Drill D2,8x13 mm DLC	Stainless Steel
	CD-GS1628	GS Drill D2,8x16 mm DLC	Stainless Steel
	CD-GS0834	GS Drill D3,4x8 mm DLC	Stainless Steel
	CD-GS1034	GS Drill D3,4x10 mm DLC	Stainless Steel
	CD-GS1134	GS Drill D3,4x11,5 mm DLC	Stainless Steel
	CD-GS1334	GS Drill D3,4x13 mm DLC	Stainless Steel
	CD-GS1634	GS Drill D3,4x16 mm DLC	Stainless Steel
	CD-GS0837	GS Drill D3,7x8 mm DLC	Stainless Steel
	CD-GS1037	GS Drill D3,7x10 mm DLC	Stainless Steel
	CD-GS1137	GS Drill D3,7x11,5 mm DLC	Stainless Steel
	CD-GS1337	GS Drill D3,7x13 mm DLC	Stainless Steel
	CD-GS1637	GS Drill D3,7x16 mm DLC	Stainless Steel
	CD-GS0841	GS Drill D4,1x8 mm DLC	Stainless Steel
	CD-GS1041	GS Drill D4,1x10 mm DLC	Stainless Steel
	CD-GS1141	GS Drill D4,1x11.5 mm DLC	Stainless Steel
	CD-GS1341	GS Drill D4,1x13 mm DLC	Stainless Steel
	CD-GS1641	GS Drill D4,1x16 mm DLC	Stainless Steel
	CD-GSFP	GS Fixation pin drill	Stainless Steel
	CD-GSST52	GS Stopper L2 Ø5	Ti6Al4V-ELI
	CD-GSTCC	GS Bone mill - Chamfer Conical	Stainless Steel

	Cat.No	Product Description	Material
	CD-GSTCH	GS Bone mill - Chamfer	Stainless Steel
	CD-GSTP	GS Tissue Punch	Stainless Steel
	CO-8040D	Digital Implant Analog Hex	Ti6Al4V-ELI
	CO-GS5090	GS Implant mount INT Hex	Stainless Steel
	CO-GS7000	GS Bone mill guide INT Hex	Ti6Al4V-ELI
	CO-GSCL	Guided Cylinder for model casting Hex	Ti6Al4V-ELI
	CO-GSFPSL	Guided sleeve for fixation pin	Ti6Al4V-ELI
	CO-GSPS	Guided Pin Screw for model casting Hex	Ti6Al4V-ELI
	CO-GSS200	GS Implant mount screw Hex	Ti6Al4V-ELI
	CO-GSSL3	GS Guided sleeve (Implant) D5 H4	Ti6Al4V-ELI
	MPN-8040D	Digital Implant Analog MPN	Ti6Al4V-ELI
	MPN-GS5090	GS Implant mount MPN	Stainless Steel
	MPN-GS7000	GS Bone mill guide MPN	Ti6Al4V-ELI
	MPN-GSCL	Guided Cylinder for model casting MPN	Ti6Al4V-ELI
	MPN-GSPS	Guided Pin Screw for model casting MPN	Ti6Al4V-ELI
	MPN-GSS200	GS Implant mount screw MPN	Ti6Al4V-ELI
	MPR-8040D	Digital Implant Analog MPR	Ti6Al4V-ELI
	MPR-GS5090	GS Implant mount MPR	Stainless Steel
	MPR-GS7000	GS Bone mill guide MPR	Ti6Al4V-ELI
	MPR-GSCL	Guided Cylinder for model casting MPR	Ti6Al4V-ELI
	MPR-GSPS	Guided Pin Screw for model casting MPR	Ti6Al4V-ELI
	MPR-GSS200	GS Implant mount screw MPR	Ti6Al4V-ELI
	MPW-8040D	Digital Implant Analog MPW	Ti6Al4V-ELI
	MU-8041D	Digital Implant Analog MU	Ti6Al4V-ELI
	CT-GS262	GS Extractor for Implant Mount	Stainless Steel
	CT-GS400	GS Motor Mount to Implant mount	Stainless Steel
	CT-GS440	GS Implant mount Extender	Stainless Steel
	CT-GSS440	GS Implant mount Short Extender	Stainless Steel
	CT-GSFP	GS Fixation pin	Ti6Al4V-ELI
	CM-0031	Square Driver One piece	Stainless Steel
	CM-0032	1 Piece adaptor	Stainless Steel
	CM-0041	Square Basic Driver Mount for OP	Stainless Steel
	CT-0219	Prosthetic Hand Driver Hex=1.25mm L7mm	Stainless Steel
	CT-0220	Prosthetic Hand Driver Hex=1.25mm L10mm	Stainless Steel

	Cat.No	Product Description	Material
	CT-0224	Prosthetic Hand Driver Hex=1.25mm L18mm	Stainless Steel
	CT-0230	Hand Ratchet	Stainless Steel
	CT-0231	Driver adaptor for hexed ratchet	Stainless Steel
	CT-0232	Driver adaptor for square ratchet	Stainless Steel
	CT-0242	Holding key hex. 2.42mm	Stainless Steel
	CT-0250	Motor Mount 1.25 Hex L22mm	Stainless Steel
	CT-0252	Motor Mount 2.42 Hex L20mm	Stainless Steel
	CT-0254	Motor Mount 2.42 Hex L28mm	Stainless Steel
	CT-0260	Drill Extension	Stainless Steel
	CT-0410	Guide pin 10mm	Stainless Steel
	CT-0413	Guide pin 13mm	Stainless Steel
	CT-0510	Direct handle driver	Stainless Steel
	CT-0510C	Direct Handle for HandPiece driver	Stainless Steel
	CT-0705	Distance indicator tool	Stainless Steel
	CT-0710	Depth Probe	Stainless Steel
	CT-0801	Ratchet with Titanium handle	Stainless Steel
	CT-0861	Torque Ratchet (40 N/cm)	Stainless Steel
	CT-0863	Torque Ratchet (100 N/cm)	Stainless Steel
	CT-D110	HandPiece Screw Driver 1.25mm L10mm	Stainless Steel
	CT-D115	HandPiece Screw Driver 1.25mm L15mm	Stainless Steel
	CT-D120	HandPiece Screw Driver 1.25mm L20mm	Stainless Steel
	CT-H207	HandPiece Implant Driver 2.42 L7mm for Hexed	Stainless Steel
	CT-H210	HandPiece Implant Driver 2.42 L10mm for Hexed	Stainless Steel
	CT-H215	HandPiece Implant Driver 2.42 L15mm for Hexed	Stainless Steel
	CT-H220	HandPiece Implant Driver 2.42 L20mm for Hexed	Stainless Steel
	CT-MU00	Multi Unit Straight abutment tighten key 10mm	Stainless Steel
	CT-MU01	Multi Unit Straight abutment tighten key 5mm	Stainless Steel
	CT-MU02	Multi Unit Straight abutment tighten key 15mm	Stainless Steel
	CT-S107	Square Driver 1.25mm L7mm	Stainless Steel
	CT-S110	Square Driver 1.25mm L10mm	Stainless Steel
	CT-S115	Square Driver 1.25mm L15mm	Stainless Steel
	CT-S207	Square Driver 2.42mm L7mm	Stainless Steel
	CT-S210	Square Driver 2.42mm L10mm	Stainless Steel
	CT-S215	Square Driver 2.42mm L15mm	Stainless Steel
	CT-S218	Square Driver 2.42mm L18mm	Stainless Steel
	CT-S262	Square abutment Extractor 26mm	Stainless Steel

Cat.No	Product Description	Material
MCT-N207	HandPiece Implant Driver L7mm for Narrow	Stainless Steel
MCT-N210	HandPiece Implant Driver L10mm for Narrow	Stainless Steel
MCT-N215	HandPiece Implant Driver L15mm for Narrow	Stainless Steel
MCT-N220	HandPiece Implant Driver L20mm for Narrow	Stainless Steel
MCT-R207	HandPiece Implant Driver L7mm for Regular	Stainless Steel
MCT-R210	HandPiece Implant Driver L10mm for Regular	Stainless Steel
MCT-R215	HandPiece Implant Driver L15mm for Regular	Stainless Steel
MCT-R220	HandPiece Implant Driver L20mm for Regular	Stainless Steel
MCT-W207	HandPiece Implant Driver L7mm for Wide	Stainless Steel
MCT-W210	HandPiece Implant Driver L10mm for Wide	Stainless Steel
MCT-W215	HandPiece Implant Driver L15mm for Wide	Stainless Steel
MCT-W220	HandPiece Implant Driver L20mm for Wide	Stainless Steel

Guided Surgical Kit

Guided Surgical Kit - CK-GS11					
Qty.	Cat. Number	Qty.	Cat. Number	Qty.	Cat. Number
1	CD-GS0822	1	CD-GS1334	1	CT-0232
1	CD-GS0828	1	CD-GS1337	1	CT-0245
1	CD-GS0834	1	CD-GS1341	1	CT-0863
1	CD-GS0837	1	CD-GS1622	1	CT-D110
1	CD-GS0841	1	CD-GS1628	1	CT-GS262
1	CD-GS1001	1	CD-GS1634	1	CT-GS400
1	CD-GS1022	1	CD-GS1637	1	CT-GS440
1	CD-GS1028	1	CD-GS1641	4	CT-GSFP
1	CD-GS1034	2	CD-GSFP	1	CT-GSS440
1	CD-GS1037	4	CD-GSST52	1	CT-H215
1	CD-GS1041	1	CD-GSTCC	1	MCT-N215
1	CD-GS1122	1	CD-GSTCH	1	MCT-R215
1	CD-GS1128	1	CD-GSTP	2	MPN-GSS090
1	CD-GS1134	1	CK-GSPS	2	MPN-GS7000
1	CD-GS1137	6	CO-GS5090	2	MPN-GSS200
1	CD-GS1141	4	CO-GS7000	4	MPR-GS5090
1	CD-GS1322	6	CO-GSS200	4	MPR-GS7000
1	CD-GS1328	1	CT-0220	4	MPR-GSS200



Surgical Kits - Internal Hex

CORKit Basic - CK-0015					
Qty.	Cat. Number	Qty.	Cat. Number	Qty.	Cat. Number
1	CT-0410	1	CD-C042	1	CD-1001
1	CT-0413	1	CT-S115	1	CD-5003
1	CT-0801	1	CT-S215	1	CD-5010
		1	CT-0260	1	CD-C038

* The CORkit basic surgical kit is used for placement of all Cortex's implants with diameters 3 mm, 3.3 mm, 3.8mm, 4.2 mm.

* For larger diameters, additional drills are required.



CORKit Prime - CK-0011					
Qty.	Cat. Number	Qty.	Cat. Number	Qty.	Cat. Number
1	CT-0252	1	CD-C050	1	CD-1001
1	CT-0254	1	CD-C060	1	CD-5003
1	CT-0260	1	CT-S107	1	CD-5010
1	CT-S262	1	CT-S115	1	CD-5137
2	CT-0410	1	CT-S207	1	CD-5140
1	CT-0801	1	CT-S215	1	CD-5148
		1	CT-0242	1	CD-C038
		1	CT-0250	1	CD-C042

* CORkit Prime surgical kit is used for placement of all Cortex's implants.



CORKit Premium - CK-0022					
Qty.	Cat. Number	Qty.	Cat. Number	Qty.	Cat. Number
1	CT-0260	1	CT-S107	1	CD-5003
1	CT-S262	1	CT-S115	1	CD-5010
2	CT-0410	1	CT-S207	1	CD-5137
1	CT-0705	1	CT-S215	1	CD-5140
1	CT-0510	1	CT-0242	1	CD-5148
1	CT-0710	1	CT-0250	1	CD-C038
1	CT-0801	1	CT-0252	1	CD-C042
1	CT-0413	1	CT-0254	1	CD-C050
1	CD-C060	1	CD-1001		

* CORkit Prime surgical kit is used for placement of all Cortex's implants.



Surgical Kits - Conical Connection

CORKit Basic - CK-0042					
Qty.	Cat. Number	Qty.	Cat. Number	Qty.	Cat. Number
1	CT-0410	1	CD-C042	1	CT-0231
1	CT-0413	1	CD-5010	1	MCT-N215
1	CD-1001	1	CT-D115	1	MCT-R215
1	CT-0260	1	CD-5003	1	MCT-W215
1	CT-0801	1	CD-C038		

* The CORkit basic surgical kit is used for placement of all Cortex's implants with diameters 3 mm, 3.3 mm, 4.2 mm.

* For larger diameters, additional drills are required.



Stoppers Kit

Prime CK-0009	
6 Stoppers for 2 mm	
Description	Cat. Number
Stopper L 6 mm for drill 2 mm	CD-0206
Stopper L 8 mm for drill 2 mm	CD-0208
Stopper L 10 mm for drill 2 mm	CD-0210
Stopper L 11.5 mm for drill 2 mm	CD-0211
Stopper L 13 mm for drill 2 mm	CD-0213
Stopper L 16 mm for drill 2 mm	CD-0216
6 Stoppers for 2.8 mm	
Stopper L 6 mm for drill 2.8 mm	CD-2806
Stopper L 8 mm for drill 2.8 mm	CD-2808
Stopper L 10 mm for drill 2.8 mm	CD-2810
Stopper L 11.5 mm for drill 2.8 mm	CD-2811
Stopper L 13 mm for drill 2.8 mm	CD-2813
Stopper L 16 mm for drill 2.8 mm	CD-2816



Premium CK-0010	
6 Stoppers for 2 mm	
Description	Cat. Number
Stopper L 6 mm for drill 2 mm	CD-0206
Stopper L 8 mm for drill 2 mm	CD-0208
Stopper L 10 mm for drill 2 mm	CD-0210
Stopper L 11.5 mm for drill 2 mm	CD-0211
Stopper L 13 mm for drill 2 mm	CD-0213
Stopper L 16 mm for drill 2 mm	CD-0216
6 Stoppers for 2.8 mm	
Stopper L 6 mm for drill 2.8 mm	CD-2806
Stopper L 8 mm for drill 2.8 mm	CD-2808
Stopper L 10 mm for drill 2.8 mm	CD-2810
Stopper L 11.5 mm for drill 2.8 mm	CD-2811
Stopper L 13 mm for drill 2.8 mm	CD-2813
Stopper L 16 mm for drill 2.8 mm	CD-2816



* Comes with two drills - pilot drill 2 mm and drill 2.8 mm.

Note

We are shaping the future of dentistry, Together.

As medical professionals, we continually seek out the next breakthrough.

A new way to practice top-quality dentistry that is, at once, innovative, personalized, user-friendly and highly effective.

*With extensive background, proven expertise far-sighted vision and pioneering digital capabilities: **We are Cortex.***

Your partner in creating a new era in dental medicine.



CORTEX[™]
The Future of Dental Implants

www.cortex-dental.com